

Product Datasheet

MYCN antibody (orb386652)

Description	Rabbit polyclonal antibody to MYCN
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human MYCN
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	50
Uniprot ID	P04198
Expiration Date	12 months from date of receipt.