

Product Datasheet

TET3 antibody (orb386661)

Description	Mouse monoclonal antibody to TET3
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human TET3
Preservatives	0.02% Thimerosal
Storage	-20°C
Note	For research use only
Clonality	Monoclonal
Purity	Caprylic Acid- Ammonium Sulfate precipitation
MW	179
Uniprot ID	043151
Expiration Date	12 months from date of receipt.