

Product Datasheet

MUC16 antibody (orb386672)

Description	Rabbit polyclonal antibody to MUC16
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human MUC16
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	2353
Uniprot ID	Q8WXI7
Expiration Date	12 months from date of receipt.