

Product Datasheet

BCL11B antibody pair (orb386675)

Description	Rabbit polyclonal antibody to BCL11B
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal [BCL11B (N) (R2122-1)] or the C-terminal [BCL11B (C) (R2122-2)]
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	96
Uniprot ID	Q9C0K0
Expiration Date	12 months from date of receipt.