

Product Datasheet

PELI1 antibody (orb386682)

Description	Rabbit polyclonal antibody to PELI1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human PELI1
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	46
Uniprot ID	Q96FA3
Expiration Date	12 months from date of receipt.