

Product Datasheet

ANKRD17 antibody (orb386684)

Description	Rabbit polyclonal antibody to ANKRD17
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human ANKRD17
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	274
Uniprot ID	075179
Expiration Date	12 months from date of receipt.