

Product Datasheet

NFKB2 antibody (orb386690)

Description	Rabbit polyclonal antibody to NFKB2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human NFKB2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	97, 52
Uniprot ID	Q00653
Expiration Date	12 months from date of receipt.