

Product Datasheet

ARX antibody (orb386702)

Description	Rabbit polyclonal antibody to ARX
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N region of human ARX
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	58
Uniprot ID	Q96QS3
Expiration Date	12 months from date of receipt.