

Product Datasheet

CHST15 antibody (orb386703)

Description	Rabbit polyclonal antibody to CHST15
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human CHST15
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	65
Uniprot ID	Q7LFX5
Expiration Date	12 months from date of receipt.