

Product Datasheet

ZNRF3 antibody (orb386706)

Description	Rabbit polyclonal antibody to ZNRF3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	C-terminal region of human ZNRF3
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	101
Uniprot ID	Q9ULT6
Expiration Date	12 months from date of receipt.