

Product Datasheet

PLEKHB2 antibody (orb386715)

Description	Rabbit polyclonal antibody to PLEKHB2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human PLEKHB2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	25
Uniprot ID	Q96CS7
Expiration Date	12 months from date of receipt.