

Product Datasheet

APOBEC3D antibody (orb386725)

Description	Rabbit polyclonal antibody to APOBEC3D
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human APOBEC3D
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	46
Uniprot ID	Q96AK3
Expiration Date	12 months from date of receipt.