

Product Datasheet

RE RE antibody (orb386729)

Description	Rabbit polyclonal antibody to RERE
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	C-terminal region of human RERE
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	172
Uniprot ID	Q9P2R6
Expiration Date	12 months from date of receipt.