

Product Datasheet

PDE3B antibody (orb386737)

Description	Rabbit polyclonal antibody to PDE3B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human PDE3B
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	124
Uniprot ID	Q13370
Expiration Date	12 months from date of receipt.