

Product Datasheet

SENP6 antibody (orb386740)

Description	Rabbit polyclonal antibody to SENP6
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human SENP6
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	126
Uniprot ID	Q9GZR1
Expiration Date	12 months from date of receipt.