

Product Datasheet

GFP antibody (orb386746)

Description	Mouse monoclonal antibody to GFP
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle terminal region of GFP from Aequorea victoria (Jellyfish)
Preservatives	0.02% Thimerosal
Storage	-20°C
Note	For research use only
Clonality	Monoclonal
Purity	Caprylic Acid- Ammonium Sulfate precipitation
MW	27
Uniprot ID	P42212
Expiration Date	12 months from date of receipt.