

Product Datasheet

MED13 antibody (orb386748)

Description	Rabbit polyclonal antibody to MED13
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human MED13
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	239
Uniprot ID	Q9UHV7
Expiration Date	12 months from date of receipt.