

Product Datasheet

BCL6 antibody (orb386751)

Description	Rabbit polyclonal antibody to BCL6
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human BCL6
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	79
Uniprot ID	P41182
Expiration Date	12 months from date of receipt.