

Product Datasheet

TET3 antibody (orb386752)

Description	Rabbit polyclonal antibody to TET3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	various regions of human or mouse TET3
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	179
Uniprot ID	O43151, Q8BG87
Expiration Date	12 months from date of receipt.