

## Product Datasheet

### NPAS4 antibody (orb386765)

|                            |                                     |
|----------------------------|-------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to NPAS4 |
| <b>Species/Host</b>        | Rabbit                              |
| <b>Reactivity</b>          | Human, Mouse                        |
| <b>Conjugation</b>         | Unconjugated                        |
| <b>Tested Applications</b> | WB                                  |
| <b>Immunogen</b>           | C-terminal region of mouse NPAS4    |
| <b>Preservatives</b>       | 0.01% NaN <sub>3</sub>              |
| <b>Storage</b>             | -20°C                               |
| <b>Note</b>                | For research use only               |
| <b>Clonality</b>           | Polyclonal                          |
| <b>Purity</b>              | Affinity purification               |
| <b>MW</b>                  | 87                                  |
| <b>Uniprot ID</b>          | <b>Q8BGD7</b>                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.     |