

Product Datasheet

TBX3 antibody (orb386767)

Description	Rabbit polyclonal antibody to TBX3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human TBX3
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	79
Uniprot ID	O15119
Expiration Date	12 months from date of receipt.