

Product Datasheet

PTPN12 antibody (orb386768)

Description	Rabbit polyclonal antibody to PTPN12
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human PTPN12
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	88
Uniprot ID	Q05209
Expiration Date	12 months from date of receipt.