

Product Datasheet

TSHR antibody (orb386777)

Description	Rabbit polyclonal antibody to TSHR
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	C-terminal region of human TSHR
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	87
Uniprot ID	P16473
Expiration Date	12 months from date of receipt.