

Product Datasheet

ZDHHC13 antibody (orb386823)

Description	Rabbit polyclonal antibody to ZDHHC13
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of mouse ZDHHC13
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	71
Uniprot ID	Q9CWU2
Expiration Date	12 months from date of receipt.