

Product Datasheet

SYK antibody (orb386832)

Description	Rabbit polyclonal antibody to SYK
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human SYK
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	72
Uniprot ID	P43405
Expiration Date	12 months from date of receipt.