

Product Datasheet

COMMD2 antibody (orb386837)

Description	Rabbit polyclonal antibody to COMMD2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human COMMD2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	23
Uniprot ID	Q86X83
Expiration Date	12 months from date of receipt.