

Product Datasheet

SLC29A3 antibody (orb386848)

Description	Rabbit polyclonal antibody to SLC29A3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SLC29A3
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	36, 52
Uniprot ID	Q9BZD2
Expiration Date	12 months from date of receipt.