

Product Datasheet

SMUG1 antibody (orb386852)

Description	Rabbit polyclonal antibody to SMUG1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human SMUG1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	31
Uniprot ID	Q53HV7
Expiration Date	12 months from date of receipt.