

Product Datasheet

RUNX1T1 antibody (orb386860)

Description	Rabbit polyclonal antibody to RUNX1T1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human RUNX1T1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	68
Uniprot ID	Q06455
Expiration Date	12 months from date of receipt.