

Product Datasheet

SEPT9 antibody (orb386893)

Description	Rabbit polyclonal antibody to SEPT9
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SEPT9
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	65
Uniprot ID	Q9UHD8
Expiration Date	12 months from date of receipt.