

Product Datasheet

WDR24 antibody (orb386902)

Description	Rabbit polyclonal antibody to WDR24
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human WDR24
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	102
Uniprot ID	Q96S15
Expiration Date	12 months from date of receipt.