

Product Datasheet

CPEB2 antibody (orb386903)

Description	Rabbit polyclonal antibody to CPEB2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human CPEB2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	65
Uniprot ID	Q7Z5Q1
Expiration Date	12 months from date of receipt.