

Product Datasheet

PANX1 antibody (orb386905)

Description	Rabbit polyclonal antibody to PANX1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human PANX1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	48
Uniprot ID	Q96RD7
Expiration Date	12 months from date of receipt.