

Product Datasheet

BAZ1A antibody (orb386914)

Description	Rabbit polyclonal antibody to BAZ1A
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human BAZ1A
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	179
Uniprot ID	Q9NRL2
Expiration Date	12 months from date of receipt.