

Product Datasheet

RECK antibody (orb386919)

Description	Rabbit polyclonal antibody to RECK
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human RECK
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	106
Uniprot ID	O95980
Expiration Date	12 months from date of receipt.