

Product Datasheet

SPO11 antibody (orb386927)

Description	Rabbit polyclonal antibody to SPO11
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	C-terminal region of human SPO11
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	45
Uniprot ID	Q9Y5K1
Expiration Date	12 months from date of receipt.