

Product Datasheet

NFKB1 antibody (orb386956)

Description	Rabbit polyclonal antibody to NFKB1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human NFKB1
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	105
Uniprot ID	P19838
Expiration Date	12 months from date of receipt.