

Product Datasheet

SRSF6 antibody (orb386990)

Description	Rabbit polyclonal antibody to SRSF6
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human SRSF6
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	40
Uniprot ID	Q13247
Expiration Date	12 months from date of receipt.