

Product Datasheet

PIK3C2A antibody (orb386999)

Description	Rabbit polyclonal antibody to PIK3C2A
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human PIK3C2A
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	191
Uniprot ID	000443
Expiration Date	12 months from date of receipt.