

Product Datasheet

CCDC6 antibody (orb387015)

Description	Rabbit polyclonal antibody to CCDC6
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human CCDC6
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	53
Uniprot ID	Q16204
Expiration Date	12 months from date of receipt.