

Product Datasheet

MAP2K2 antibody (orb387028)

Description	Rabbit polyclonal antibody to MAP2K2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human MAP2K2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	44
Uniprot ID	P36507
Expiration Date	12 months from date of receipt.