

Product Datasheet

TLR5 antibody (orb387033)

Description	Rabbit polyclonal antibody to TLR5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human TLR5
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	98
Uniprot ID	O60602
Expiration Date	12 months from date of receipt.