

Product Datasheet

Tubulin A antibody (orb387040)

Description	Rabbit polyclonal antibody to Tubulin A
Species/Host	Rabbit
Reactivity	Human, Mouse, Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of Arabidopsis Tubulin A1 (At1g64740)
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	49
Uniprot ID	P11139
Expiration Date	12 months from date of receipt.