

Product Datasheet

TFG antibody (orb387043)

Description	Rabbit polyclonal antibody to TFG
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human TFG
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	43
Uniprot ID	Q92734
Expiration Date	12 months from date of receipt.