

Product Datasheet

TBC1D14 antibody (Biotin) (orb38706)

Description	TBC1D14 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human TBC1 domain family member 14 protein (136-435AA)
Target	TBC1D14
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9P2M4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.