

Product Datasheet

CHD8 antibody (orb387078)

Description	Rabbit polyclonal antibody to CHD8
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human CHD8
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	291
Uniprot ID	Q9HCK8
Expiration Date	12 months from date of receipt.