


Product Datasheet

ZNF804A antibody (orb387083)

Description	Rabbit polyclonal antibody to ZNF804A
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human ZNF804A
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	137
Uniprot ID	Q7Z570
Dilution Range	1:1,000-1:3,000
Expiration Date	12 months from date of receipt.



Western blot analysis of A) mouse thymus...