

Product Datasheet

TLR7 antibody (orb387087)

Description	Rabbit polyclonal antibody to TLR7
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human TLR7
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	121
Uniprot ID	Q9NYK1
Expiration Date	12 months from date of receipt.