

Product Datasheet

TICAM2 antibody (orb387091)

Description	Rabbit polyclonal antibody to TICAM2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human TICAM2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	27, 46
Uniprot ID	Q86XR7
Expiration Date	12 months from date of receipt.