

## Product Datasheet

### KIRREL3 antibody (orb387102)

|                            |                                       |
|----------------------------|---------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to KIRREL3 |
| <b>Species/Host</b>        | Rabbit                                |
| <b>Reactivity</b>          | Human                                 |
| <b>Conjugation</b>         | Unconjugated                          |
| <b>Tested Applications</b> | WB                                    |
| <b>Immunogen</b>           | N-terminal region of human KIRREL3    |
| <b>Preservatives</b>       | 0.01% NaN <sub>3</sub>                |
| <b>Storage</b>             | -20°C                                 |
| <b>Note</b>                | For research use only                 |
| <b>Clonality</b>           | Polyclonal                            |
| <b>Purity</b>              | Affinity purification                 |
| <b>MW</b>                  | 85                                    |
| <b>Uniprot ID</b>          | <b>Q8IZU9</b>                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.       |